










<div><div><div><div><div></div></div></div><div>LABORATORIES</div></div><div>TEOXANE</div><div>GENEVA</div><div>CONFIDENT IN YOUR SKIN. CONFIDENT IN YOUR LIFE™</div></div>	TEOSYAL® DERMAL FILLERS RANGE									
	THE ENHANCERS		THE ESSENTIALS				THE VOLUMIZERS		THE SPECIFICS	
	Meso	Redensity [I]	First Lines	Touch Up	Global Action	Deep Lines	Ultra Deep	Ultimate	Kiss	Redensity [II]
HA Concentration (mg/ml)	15	15	20	25		25	25	22	25	15
Lidocaine concentration for Teosyal® PureSense Line	no lidocaine	0,3% <sup>(2)</sup>	0,3%	no lidocaine	0,3%	0,3%	0,3%	0,3%	0,3%	0,3% <sup>(2)</sup>
Composition	Non cross-linked HA	Non cross-linked HA + Dermo-Restructuring Complex <sup>(1)</sup>	Cross-linked HA	Cross-linked HA		Cross-linked HA	Cross-linked HA	Cross-linked HA	Cross-linked HA	Mix of non cross-linked and cross-linked HA + Dermo-Restructuring Complex <sup>(1)</sup>
Cross-linking level	○ ○ ○ ○ ○ ○	○ ○ ○ ○ ○ ○	● ○ ○ ○ ○ ○	● ● ● ○ ○ ○		● ● ● ● ○ ○	● ● ● ● ● ○	● ● ● ● ● ●	● ● ● ● ○ ○	● ○ ○ ○ ○ ○
Treated Areas and Indications	Hydration of the face, neck and décolleté	Redensification and hydration of the face, neck and décolleté	Correction of fine lines and wrinkles in thin skin (crow's feet)	Ideal for touch up as complement of a first Teosyal® injection	Filling of moderately pronounced wrinkles	Filling of highly pronounced wrinkles	Volume creation for targeted areas (cheekbones, chin) and correction of deep wrinkles in thick skin	Volume recovery for wide areas (face oval, cheeks, temples)	Harmonization of the lip contour and volume	Correction of the eye circle, tear trough, palpebromalar groove, minimal palpebral or malar bags
Injection Depth	Superficial dermis 	Superficial to mid dermis 	Superficial dermis 	Mid dermis 		Deep Dermis 	Subcutaneous and/or pre-periosteum 	Subcutaneous or pre-periosteum 	Labial mucosa 	Supra-periosteal, deep sub-orbicular 
Blister pack	2 x 1ml	1 x 3ml 2 x 1ml	2 x 0,7ml	2 x 0.5ml	2 x 1ml	2 x 1ml	2 x 1.2ml	1 x 3ml 2 x 1ml	2 x 1ml	2 x 1ml
Needle	30G <sup>1/2"</sup>	30G <sup>1/2"</sup>	30G <sup>1/2"</sup>	30G <sup>1/2"</sup>		27G <sup>1/2"</sup>	25G <sup>1"</sup>	23G <sup>1"</sup> 27G <sup>1/2"</sup> <sup>(3)</sup>	27G <sup>1/2"</sup>	30G <sup>1/2"</sup>
Estimated duration <sup>(4)</sup>	2 - 3 months	Protocol: 3 sessions every 3 weeks, twice a year	6 months	Up to 9 months		9 months	12 months	Up to 18 months	Up to 9 months	9 - 12 months

(1) Dermo-Restructuring Complex composition: 8 amino acids: Glycine, Lysine, Threonine, Proline, Isoleucine, Leucine, Valine, Arginine, 3 antioxydants: Glutathione, N-acetyl-L-cysteine, α-Lipoic acid, 2 minerals: Zinc, Copper, 1 vitamin: B6

(2) Redensity I and Redensity II references are only available in the Teosyal® PureSense range


(3) for Ultimate 2 x 1 ml

(4) Average duration depends on several factors: patient's skin type, severity of the wrinkle to be corrected, type of injection and volume injected

# TEOSYAL


EFFECTIVE • SAFE • LONG-LASTING

The Enhancers ■ The Essentials ■ The Volumizers ■ The Specifics



Eyes

The Specifics



Lips

  
LABORATORIES  
**TEOXANE**  
GENEVA

CONFIDENT IN YOUR SKIN, CONFIDENT IN YOUR LIFE™